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Instruction manual lightweight wheelchair

Artikelnr.: KREHA-2500-43

GB Vorwort

Hello Customers!

Thank you for your love of Yuyue and for purchasing a Yuyue Electric Wheelchair.

The electric wheelchair carefully developed by Yuyue is light and convenient, efficient and energy-saving. It is flexible and safe to operate, and is well received by the majority of users. Before use, please carefully read these instructions so that you can better understand various functions of the electric wheelchair for better control. You should carry out the maintenance and service as required to ensure that the wheelchair is in good condition. For any questions, please contact the dealer or manufacturer and we will serve you wholeheartedly!

Safety Overview

- ⓘ Do not operate your electric wheelchair before reading and fully understanding these instructions.
- ⓘ Do not operate your electric wheelchair before the assembly and inspection are completed.
- ⓘ It is recommended that persons with unsound minds, slow responses and operational difficulties should not use the electric wheelchair.
- ⓘ Do not disassemble or modify the electric wheelchair or use any replacement parts not manufactured by the company.
- ⓘ Do not get into or out of the electric wheelchair when the controller is

powered on or the electric wheelchair is in the manual mode and is not fixed by the pusher.

- ⓘ Do not use the electric wheelchair when the anti-roll wheel is not open or has broken down.
- ⓘ Do not tilt or lift the electric wheelchair to one side when it is being operated normally.
- ⓘ Do not stand on the pedal to prevent the electric wheelchair from rolling over.
- ⓘ Do not turn or steer on an incline. To avoid accidents, the electric wheelchair should not be operated by two persons simultaneously.
- ⓘ The moving electric wheelchair should be decelerated to below 2km/h before turning.
- ⓘ Drive carefully in the lowest gear when going downhill.
- ⓘ Please check whether the wheel connections are secure and reliable.
Pull the controller joystick gently and do not pull back and forth on it rapidly.
- ⓘ The controller is the core part of the electric wheelchair. Do not park the electric wheelchair in the open air for a long time. The electric wheelchair should be parked indoors when it rains, to avoid moisture.
- ⓘ Before powering on the controller, confirm that the "Manual/Electric" clutch of the left and right motors or the handle is in the "Electric" gear.

- ① Do not switch the "Manual/Electric" clutch handle of the motor to the "Manual" position while the wheelchair is running; if it is necessary to use radio communication devices such as mobile phones or laptops while in the electric state, please turn off the controller power of the electric wheelchair.
- ① The electric wheelchair is suitable for flat ground and low inclines. Avoid driving on pavement with a slope greater than 6 degrees and driving over obstacles more than 4cm high.
- ① It is strictly forbidden to cross places with horizontal gaps, such as sewers.
- ① If the electric wheelchair is not used for a long time, turn off the power switch on the battery box.
- ① This company's electric wheelchair is suitable for using outdoors, but can only be used in neighborhoods.
- ① This electric wheelchair shall not be operated on the road, otherwise it is easy to cause traffic accidents.



3. Product Characteristics

(1) Scope of application

This company's electric wheelchair is suitable for the disabled with difficulties walking and the elderly and infirm.

(2) Product components

The product is composed of a frame, controller, motor, battery, pedal, handrail, front wheels and rear wheels.

GB

(3) Structural characteristics

- Motor: energy-saving and efficient.
- Frame: easy to assemble and disassemble, foldable
- Intelligent controllers: power button, power display, universal joystick, horn
- Electromagnetic braking system: safe and reliable
- Anti-roll device
- Battery: fully-sealed, maintenance-free
- Two modes of operation: automatic electric drive mode and powerboost mode

(4) Technical parameters

- Product type: outdoor model
- Max. speed: ≤ 6 km/h
- Ambient operating temperature: $-25^{\circ}\text{C} \sim +50^{\circ}\text{C}$
- Power consumption: $\leq 2,5$ kwh/100km
- Load capacity: ≤ 100 kg
- Braking performance on horizontal plane: $\leq 1,5$ m
- Endurance mileage: ca. 20km
- Max. safety slope braking: $\leq 3,6$ m (6°)
- Battery: lithium battery DC 24Vx18Ah; lead acid battery DC 24Vx21Ah
- Obstacle clearance height: 40mm
- Groove clearance width: 10 cm
- Gradeability: $\leq 6^{\circ}$
- Hill-holding performance: 9°
- Static stability: $\leq 9^{\circ}$
- Dynamic stability: $\geq 6^{\circ}$
- Min. turning radius: 1,2m
- Noise: < 65 dB

The above parameters will vary depending on the weight of the occupant, road conditions and battery usage

- Normal working conditions:

Ambient temperature range:

$-25^{\circ}\text{C} \sim +50^{\circ}\text{C}$

Relative humidity range:

25%~95%

Atmospheric pressure range:

86kPa~106kPa

- Internal power supply: DC24V \pm 5V
- Internal power supply equipment
- Electrical requirements: type B application part
- Mode of operation: continuous operation
- Fluid ingress protection degree: IPX4
- Motor power: about 130W for both left and right motors
- This equipment cannot be used in flammable anesthetic gas mixed with air or flammable anesthetic gas mixed with oxygen or nitrous oxide.

Daily Maintenance and Service

Hinweis:

Before servicing, press the power switch on the battery box to the "O" state.

- It is strongly recommended to have specialized personnel adjust and replace worn parts or to contact the manufacturer.

| Inspektion intervall | Daily | Weekly | Monthly | Every 3 months | Every 3 months |
|--------------------------|-------|--------|---------|-------------------|-------------------|
| Batterie | • | | | | |
| Tire pressure | | • | | | |
| Wire | | • | | | |
| Hand brake | | | • | | |
| Frame | | | | • | |
| Controller | | | | • | |
| Motor | | | | • | |
| Connector | | | | • | |
| Seat back cushion | | | | | • |
| Tire | | | | | • |
| Electromagnetic brake | | | | | • |

GB

Worn parts are replaced as follows
(if it is difficult to replace, please contact the manufacturer promptly for a replacement):

1. Replacement of front wheel:

Loosen the screws with a wrench, remove the front wheel, install the new front wheel, lock the screws, adjust the screw tightness and make sure the front wheel spins smoothly.

2. Replacement of rear wheel:

have specialized personnel replace it or contact the manufacturer.

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Loosen the screws with a screwdriver, remove the seat (back) cushion, install the new seat (back) cushion and lock the screws with a screwdriver.

4. Replacement of handrail:

Loosen the screws with a wrench, remove and handrail, fasten the new handrail and lock the screws with a wrench.

• Battery:

Check the remaining battery capacity. If the battery reaches the end of its service life and should be replaced, please

contact the supplier or the manufacturer or purchase a battery of the same specification locally.

• Tire pressure:

It is recommended to inflate the tires of 200 X 45 - 110 and 310 X 50 - 203 to 260 kpa (maximum 325 kpa) and inflate the tires of 22 X 1.75 to 280 kpa (maximum 345kpa), which can be adjusted depending on weight and air temperature variation; the tires will leak gradually if they are stored for a long time or are not used, this is normal; please carefully read the following instructions to operate the wheelchair: a. when the tire pressure is not sufficient, (1) press the tire evenly with your hands, so that the tire and rim can fit uniformly; (2) inflate the tire to the appropriate state; b. when the tire is deflated, (1) inflate the tire to about 30% of the full air volume, press the tire evenly with your hands, so that the tire and rim can fit uniformly; (2) inflate the tire to the appropriate state.

• Wire:

Check the electrical parts and connecting wires for any damage. If any, please contact the supplier or have specialized personnel repair it. Do not repair it by yourself.

• Hand brake:

The brake prevents the wheels from moving after the wheelchair is parked and cannot be used during driving; check whether the brake has broken down. If it has, it can be adjusted and resumed.

• Frame:

The frame surface coating should be wiped with a soft cloth and should be kept clean; lubricant is prohibited. If the frame cracks, please contact the supplier

· **Service of controller:**

Clean the controller and joystick with a cloth dampened with a neutral diluted cleaner carefully. Never use abrasive materials or alcoholic detergents for cleaning. Protect the controller from damage during transportation of the wheelchair.

· **Motor:**

Check if there is any increase in oil leakage or noise. If any, please contact the supplier or manufacturer. Maintenance of connector reliability: check that the screws and nuts on the body are tightened frequently and promptly handle the problems (if any) to ensure driving safety.

· **Seat back cushion:**

Wash the seat cover and backrest with warm water and diluted soapy water. Avoid keeping the wheelchair in a damp place.

· **Electromagnetic brake:**

Service brake. Let the wheelchair run straight at the maximum speed on a flat asphalt pavement. Release the controller joystick to automatically return it to the original position and measure the distance from the time of releasing the joystick to stopping. If the distance is larger than it was originally, the braking effect is reduced; if the distance is more than 1.5m, contact the supplier or manufacturer for repair

· **Maintenance and service of battery:**

1. Pay attention to the power indicator on the controller panel: If the display shows no green lights, you should charge the battery soon. If only the red light is on and no green or yellow lights, the charge level is extremely low. Please recharge the battery immediately to avoid low battery voltage, that affects the battery life
2. The battery is marked with obvious positive and negative electrode signs and a reliable connector is provided to ensure

the normal connection of the circuit.

Non-specialized persons should not connect the circuit.

3. The battery is maintenance-free and replenishment of supplemental fluid is not required daily. During charging, the battery temperature will rise, but should not exceed 45°C; if the temperature exceeds 45°C, stop charging until the temperature drops below 35°C. If the wheelchair is parked for a long time, recharge the battery at least once a month.

4. The battery has its service life. If the trip mileage is significantly different from the nominal mileage after long time of normal operation, please replace the battery.

5. Do not use the battery at temperature of $\geq 50^{\circ}\text{C}$ oder $\leq -20^{\circ}\text{C}$.

6. Keep the battery clean and dry. Do not hit the battery with hard objects. Keep the battery properly and keep it in a place out of reach of children.

7. The power switch on the battery box cuts off the battery power and reduces the natural power consumption of the battery. Turn off the power switch on the battery box when the wheelchair is not being used for a long time.

8. "Full capacity": develop the habit of maintaining the battery at full capacity and charge the battery in a timely manner according to usage, so that the battery is at "full capacity" long-term.

9. Replace battery: Press the power switch on the battery box to the "O" position for battery replacement;

(1) If lead-acid batteries is used, left and right battery boxes are provided and each battery box contains a lead-acid battery: red on the battery indicates the positive electrode and

black indicates the negative electrode. The battery shall be replaced as follows:

- a. Loosen the 6 self-tapping screws from inside of the electric wheelchair (i.e. back of the battery box) and open the back cover, see Figure 19. The methods are same for left side and right side;


Abb 19

- b. Take the battery out of the battery box, remove the positive and negative leads with a wrench, install the new battery and affix the positive and negative leads. The remaining wiring harness shall remain unchanged;

- c. Pay attention to distinguish between the positive and negative electrodes on both sides of the wiring harness. For the right hand battery connecting wires, connect the red wire with the positive electrode (+) and connect the brown wire with the negative electrode (-); for the battery on the left, connect the red wire with the positive electrode (+) and connect the black wire with the negative electrode (-)

- (2) If a lithium battery is used, only one battery box is provided at the right and the battery box only has a lithium battery: red on the battery indicates the positive electrode and black indicates the negative electrode. The battery shall be replaced as follows:

- a. Loosen the 6 self-tapping screws from inside of the electric wheelchair (i.e. back of the battery box) and open the back cover (the illustration is same as that of the lead-acid battery);
- b. Take the battery out of the battery box, separate the connector by hand, install the new battery and replace the connector. The remaining wiring harness shall remain unchanged.
- Waste should be disposed of in accordance with the relevant national environmental protection regulations.

Transportation and Storage

(1) Transportation

- Handle the product gently during transportation. Do not throw, turn over or apply great pressure to it.
- Folding steps: (see Figures 20, 21 and 22)


Abb 20

Abb 21

Abb 22

- Before boarding a plane, run out the battery capacity until the battery capacity is low (up to 1/3 of total battery capacity) and the red or red-yellow indicator lights are lit or the voltage is within 23V (lithium battery; within 26V).
- Turn off the power switch on the battery box during transportation

(2) Storage

This product should be placed in dry and ventilated place and should not be stored in a place subjected to high temperatures and rapid temperature changes; this product should be separated from acid, alkali and other chemical corrosive items.

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Fault Diagnosis and Handling

- If the power signal is not available when the power button on the controller is pressed, confirm whether the power switch on the battery box is at the "I" position; if it is at the "O" position, press the power switch on the battery box to the "I" position; then, if the power signal is still not available for the power button on the controller, the current may be too large and the overload protector has cut off the

Abb 23

power automatically. After the problem is eliminated, press the button to reset manually (see Figure 23).

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GB

- In case of failure, the controller will sound an alarm and blink, and the corresponding faults can be handled according to the times of the alarm sounds. The diagnostic number displayed by the product's built-in diagnostic function can reflect the nature of the abnormal condition. These abnormal conditions can still be detected without other service tools. Acoustic signal: the first 2 high-frequency sounds are the guide sounds and the subsequent sounds are the times of the alarm sounds in cycles.

Diagnostic description and recommended solutions:

Low voltage:

If the accumulator voltage is low, use it after charging; if the accumulator is damaged, replace it; the battery cannot be charged.

Right or left motor failure:

Check the right motor for any loose connecting devices or motor wires.

Right or left brake:

Check the right or left brake for any loose connecting devices or motor wires. Check the brake switch for any damage or poor switch contact.

Overcurrent:

Check the brake and motor transmission for any seizure. If the current detected by the ampere meter is not large, the controller may have broken down.

Rocker:

The rocker is not reset or the rocker wire is broken or the connector has become loose.

Controller failure:

Consult the service manufacturer.

| | | |
|---|--|---|
| 1 | The power indicator does not light when the controller power button is pressed | <p>The connection between the battery and the controller is not correct. Re-connect it after inspection.</p> <p>The battery voltage is too low. If the indicator still does not light up after the battery is charged, the battery may have reached the end of its service life. Please replace the battery.</p> <p>If the controller breaks down, contact the dealer, manufacturer or after-sales service department.</p> |
| 2 | Low battery voltage | <p>The connector contact between the battery and the controller is not reliable. Please insert it again.</p> <p>The contact resistance between the battery connection is large; if the contact resistance is not large, the contact surface may be oxidized or loosened. Remove the oxide layer or install the connector properly.</p> <p>If the battery voltage is too low, it indicates that the battery has reached the end of its service life. Please replace the battery.</p> |
| 3 | High battery voltage | <p>The battery's charging voltage is too large. The voltage should not exceed 29.4V after charging is completed.</p> |

GB

4

The motor
does not work

If the motor breaks down, contact the dealer, manufacturer or after-sales service department.

The motor connector is not connected reliably. Please connect it again.

5

Motor brake
failure

The motor connector is not connected reliably. Please connect it again.

The electric brake coil is damaged.

6

No charging
indication

The connector of the motor and the controller is not connected reliably. Please connect it reliably again.

The battery has reached the end of its service life or the charger is damaged. Please replace the battery or charger.

7

Trip distance is
shortened after
charging

The battery is not charged fully. Please charge the battery again.

The battery is approaching the end of its service life. Please replace the battery.

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KAPHINGST ONLINE AG
Proteno GmbH