REHASHOP Elektromobil Siam

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
D

Thank you very much for your confidence in our products.

Please read the following information and instructions carefully before using your new REHASHOP electric scooter Siam.

For a safe handling of your new scooter, please note the instructions for use and safety in particular.

If you have any questions, please contact our customer service team. You can reach them free of charge at the number **0800 888 90 80**.

We wish you all the best with your new product!

Your rehashop team

GB

Table of Content

Page

1. Import precautions	25
2. Safety information on electromagnetic interference (EMI).....	26
GB 3. Safety warning and instruction labels	28
4. Identification of parts	29
5. Operating your scooter	30
6. Charging the batteries	33
7. Disassembling your scooter	34
8. Caution: Obstacle climbing	35
9. Care and maintenance	36
10. Other information	38
11. Error code	40
12. Specification	42
13. Spare parts list REHASHOP electro scooter Siam	43

Manufacturer:

CHIEN TI ENTERPRISE CO., LTD.

No. 13, Lane 227, Fu Ying Rd., Hsin Chuang District,
New Taipei City, Taiwan

Tel : +886-2-2903-2987

Fax : +886-2-2903-8807

E-Mail : sales@chienti.com.tw

http : //www.chienti.com.tw

Information of European Representative :

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

Symbols used in this manual:


 **Symbol for warnings**

Please follow the instructions on the side of this symbol. Ignoring them may result in physical injury or damage to the Scooter and your surroundings.

GB

1. Important Precautions

- Only one person at a time could ride the REHASHOP Scooter Siam.
- Maximum load is 136 kg / 300 lbs pounds.
- Turn the key off before getting on or off.
- Always drive carefully and be aware of others using the same area.
- Always use pedestrian crossings wherever possible. Take extreme care when crossing roads.
- Do not drive on slope exceeding 8 degree and take extreme care when turning on slope.
- Do not use full power when turning to sharp corner.
- Take great care and drive in low speed when backing up, riding downhill or on uneven surface and climbing curb.
- Please use the lowest speed when driving through the descending road or uneven terrain. If speed is too fast, leave your hand off the handle bar, let the scooter stop. Make sure safety and start again.
- The weight capacity limit at different ramp degree (please refer to following picture).

*The grade Climbing degree will be affected by weight capacity, driving speed, an ramp degree, and scooter parameter.

*To prevent any danger from motor defected; please avoid to drive on long ramp or any uneven terrain.

- A slow speed must always be used when ascending, descending or traversing aslope or incline and also on uneven terrain, ramps and soft or loose surfaces, such as gravel or grass.
- To prevent any danger, do not turn around at high speed on ascending, descending ramp.
- Scooter may not operate well in high humidity.

- Do not leave the powered scooter in a rain storm of any kind.
- Direct exposure to rain or dampness will cause the scooter to malfunction electrically and mechanically; may cause the powered scooter to prematurely rust.
- Never put scooter in neutral when staying on slopes.
- Follow traffic laws when riding outside.
- When scooter on moving transport vehicles, do not sit or stay on scooter.

GB

2. Safety Information on electromagnetic interference (EMI)

It is very important that you read this information regarding the possible effects of electromagnetic interference on your motorized scooter.

Powered wheelchairs and motorized scooters may be susceptible to electromagnetic interference (EMI), which is interfering electromagnetic energy (EM) emitted from sources such as radio stations, TV stations, amateur radio (HAM) transmitters, two-way radios, and cellular phones. The interference (from radio wave sources) can cause the motorized scooter to release its brakes, move by itself, or move in unintended directions. It can also permanently damage the motorized scooter control system. The intensity of the interfering EM energy can be measured in volts per meter (V/m). Each motorized scooter can resist EMI up to certain intensity. This is called its „immunity level.“ The higher the immunity level, the greater the protection. At this time, current technology is capable of achieving at least a 20 V/m immunity level, which would provide useful protection from the more common sources of radiated EMI. This immunity level of this motorized scooter model is 20 V/m.

There are a number of sources of relatively intense electromagnetic fields in the everyday environment. Some of these sources are obvious and easy to avoid. Others are not apparent and exposure is unavoidable. However, we believe that by following the warnings listed below, your risk to EMI will be minimized.

The Sources of radiated EMI can be broadly classified into three types:

1. Hand-held portable transceivers (transmitters-receivers) with the antenna mounted directly on the transmitting unit. Examples include: citizens band (CB) radios, „walkie talkie,“ security, fire, and police transceivers, cellular telephones, and other personal communication devices;

Some cellular telephones and similar devices transmit signals while they are ON, even when not being used.

2. Medium-range mobile transceivers, such as those used in police cars, fire trucks, ambulance and taxis. These usually have the antenna mounted on the outside of the vehicle; and
3. Long-range transmitters and transceivers such as commercial broadcast transmitters (radio and TV broadcast antenna towers) and amateur (HAM) radios.

Other types of hand-held devices, such as cordless phones, laptop computers, AM/FM radios, TV sets, CD players, cassette players, small appliances, such as electric shavers and hair dryers, so far as we know, are not likely to cause EMI problems to your motorized scooter.

Motorized Scooter Electromagnetic Interference:

Because EM energy rapidly becomes more intense as one moves closer to the transmitting antenna (source), the EM fields from hand-held radio wave sources (transceivers) are of special concern. It is possible to unintentionally bring high levels of EM energy very close to the motorized scooter control system while using these devices. This can affect motorized scooter movement and braking. Therefore, the warnings listed below are recommended to prevent possible interference with the control system of the motorized scooter.

Warnings:

Electromagnetic interference (EMI) from sources such as radio and TV stations, amateur radio (HAM) transmitters, two-way radios and cellular phones can affect motorized scooters. Following the warnings listed below should reduce the chance of unintended brake release or motorized scooter movement which could result in serious injury.

1. Do not operate hand-held transceivers (transmitters-receivers), such as citizens band (CB) radios, or turn ON personal communication devices, such as cellular phones, while the motorized scooter is turned ON;
2. Be aware of nearby transmitters, such as radio or TV stations, and try to avoid coming close to them;
3. If unintended movement or brake release occurs, turn the motorized scooter OFF as soon as it is safe;
4. Be aware that adding accessories or components, or modifying the motorized scooter, may make it more susceptible to EMI;



There is no easy way to evaluate their effect on the overall immunity of the motorized scooter.

5. *Report all incidents of unintended movement or brake release to the distributor listed on the inside front cover of this manual. Note whether there is a source of EMI nearby.*

3. Safety warnings and instruction labels :

GB

	<p>Warning Sticker</p> <ol style="list-style-type: none"> 1. Please read the instruction booklet carefully before using your scooter. 2. Do not drive the scooter on slippery surfaces or on slopes over 8 degrees limit. 3. Do not drive on highway, crowded roads, or unfamiliar areas. 4. Do not turn at high speed in either forward or reverse. 5. Do no wash with water or leave scooter in humid environment since water can damage the electronic parts. 		
	<p>Power warning label</p>		<p>Information Sticker states the date of manufacture and serial number of the scooter.</p>
	<p>Wiring diagram label</p>		<p>This product has been supplied from an environmentally aware manufacturer.</p> <p>It may contain substances that could be harmful to the environment.</p>
	<p>Instruction of battery lock nut label</p>		<p>Recycling must be carried out in accordance with the respective national legal provisions.</p>
	<p>Warning label: "Keep hands away"</p>		<p>N-D Lever Adjustment label which instructs freewheel mode operation.</p>

4. Identification of Parts:

Before attempting to drive this scooter on your own, it is important that you familiarize yourself with the controls, and how to operate them.

Figure 1: Front View

- A:** Easy-To-Operate Control Panel & Delta Handlebar
- B:** Flip-up Armrests
- C:** Padded, Foldable Seat
- D:** Front Basket
- E:** Tiller Angle Adjustment
- F:** Seat Rotation Lever
- G:** Release Lever
- H:** Reflector
- I:** Charging Socket
- J:** Solid Tires

Figure 2: Top Control Panel

- 1:** Horn
- 2:** Self-Diagnostic Warning Light
- 3:** Battery Gauge
- 4:** Speed Dial Knob
- 5:** Throttle
- 6:** Main Key Switch

6

Figure 3: Back View


- 7:** Reflectors
- 8:** Anti-Tip Wheels
- 9:** Free-Wheeling Lever
- 10:** Armrest Adjustment Knob

5. Operating your Scooter :

Before operation, please remove the Lock Nut (A) and press button of Circuit Breaker (B) to activate the system.

GB

Function of Parts:

Main Key Switch (C)

1. Turn the Key to the right.
Turn the Scooter on.
2. Turn the Key to the left.
Turn the Scooter off.

Always ensure that the scooter is switched off before getting on or off the scooter and before removing any items of the scooter.

Turning the scooter OFF whilst driving will bring the scooter to an abrupt stop.

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The rabbit icon indicates fast or high speed. The turtle icon indicates slow or low speed.

By turning this Knob (D) towards chosen icon you can control overall speed of the scooter.

Before driving the scooter, set the speed to low speed by turning the knob towards turtle icon.


Drive in high speed when encounter an up slope, and drive in low speed when encounter a down slope.

Horn Button (E)

Press Horn Button (E) once to sound warning tone when necessary.

Battery Gauge (F)

The LEDs represents an approximation of battery capacity. If the charge is full, all LEDs are lit. When only red and orange LEDs are lit, the batteries are consuming and should be charged soon. When only red LEDs are lit, charge the batteries immediately.

Self Diagnostic Warning Light (G)

Flashing light indicates there is a problem within scooter. See Page 20 for more Information.

Throttle (H) (H1)

To move forward by pulling the right side of Throttle (H). To move backward by pulling the left side of Throttle (H1). (The movement directions can be reversed by local dealers if required.) Releasing the throttle to engage automatic brake. These are also your accelerator, the further you pull it, the faster you go. (Subject to the position of the Rabbit/Turtle control).

Releasing the throttle engages the automatic brake, but will taxi for a short buffering distance, please keep safe distance when stopping to prevent any danger.

Tiller Angle Adjustment:

1. Turn the Cap (I) outward upward to disengage the pin.
2. Simultaneously, adjust the tiller fore-and-aft to the most comfortable angle.
Release Cap (I) and ensure the pin is fully engaged to lock the tiller in position.

Seat rotation Adjustment:

1. Lift the Lever (J) upward to disengage the pin.
2. Simultaneously, rotate Seat (K) to the most comfortable angle. To lock the seat in position by Releasing Lever (J) and ensure the pin is fully engaged.

GB

Circuit Breaker:

If scooter's circuit system malfunctions or over loaded, the circuit breaker will trip and automatically shut down the power to ensure driver's safety. After shut down the power, press button of Circuit Breaker (B) to reactive the circuit system.

Free-Wheeling Lever:

- The scooter can be pushed, when the Free-Wheeling Lever (L) is disengaged.
- To engage the Free-Wheeling Lever (L) - Pull up the free-wheeling lever.
- To disengage the Free-Wheeling Lever (L) - Push down the free-wheeling lever.

Basket assembling and disassembling:

1. Place the Front Basket (M) to the bracket on the tiller.
2. Remove the Front Basket (M) and place it aside.

Armrest assembling and angle adjustment:

1. Adjust armrest to most comfortable length, then lock the Knob (N).
2. Adjust Screw (O) height to preferable armrest's angle.

Seat height adjustment:

1. Remove the seat, then screw, nut and washer from seat post.
2. Adjust seat post to desired height, then attached tightly with screw, nut and washers.
3. Assemble the seat back to its original position.
4. Seat height adjustability HS-295 485mm / 510mm / 535mm.

6. Charging the Batteries:

Your REHASHOP scooter is equipped with two sealed, maintenance free 12V. 12Ah. Rechargeable batteries and one 2A On-board charger. Batteries must be charged before using scooter for first time and should be recharged after each day's use. Be sure power switch is in OFF position and free-wheeling lever is not in freewheel mode.

The batteries should be charged for 24 hours before first use to ensure the longest possible life for your scooter.

The scooter should be charged at temperatures between 10°C and 30°C.

For safety reasons, only sealed non-spillable batteries that meet DOT CFR 173.159(d), IATA Packing Instructions, and IATA Provision A67 shall be installed in the scooters. If you need new batteries, contact your dealer. Only 12V. 12Ah Valve-regulated Lead Acid batteries should be used.

1. Insert battery charger cord into Charging Socket (P) on battery pack.
2. Plug other end of Power Cable (Q) into a standard electrical wall outlet.
3. The charging indicator will normally be red or yellow at this point.
4. Allow batteries to charge until charging indicator turns green.
5. After indicator turns green, unplug battery charger from scooter and wall outlet.
6. Do not charging simultaneously with battery box socket and optional charging socket, this will cause serious risk.
7. If at any time battery charger light flashes green over 40 minutes, this indicates abnormal charging occurred.

You should check the following:

- Charger plug is correctly positioned
 - Scooter is turned off
- s.o. are the problems, contact your local authorized dealer.

The time needed to recharge will vary depending on the depletion of the batteries. Charging for longer than necessary will not harm the batteries. They can not be overcharged.

Keep in Mind these Rules:

- Fully charge batteries at least once a month, or more if you use scooter regularly. Charge after each trip exceeding 3 kilometers.
- If storing your scooter for some time (1 month or more), make sure that batteries are fully charged, and on returning, charge them again before using scooter.
- Batteries will only give maximum performance after scooter has been used, and batteries have been recharged up to 10 times. A bit like breaking in a new car.

GB

Please be aware that the travelling range of your mobility scooter is impacted by how fast the batteries are discharged. This will depend on many circumstances, such as ambient temperature, condition of the surface of the road, tyre pressure, weight of the driver, driving environment (inclines etc.) and utilization of your lighting system if fitted. We recommend that you test your local ride with a family member to ensure a safe journey.

7. Disassembling your Scooter:

I. **Seat disassembling:**

Remove seat by lifting Seat Rotation Lever (J). Then rotate Seat (K) and lift up away from scooter.

II. **Battery Pack disassembling:**

As indication, pull Battery Pack Handle (R) to remove Battery Box (S) from scooter. Caution: Battery Box is heavy. When lifting, please use correct lifting posture to avoid injury. Ask for assistance if necessary.

When assembling Battery Pack (T), make sure to aim for the battery terminal connection for proper electricity conductivity.

III. **Front & Rear Frame disassembling:**

Hold front & rear frame with Release Lever (T) to disassemble front and rear frame.

Tiller fold method:

Turn the Cap (I) to fold down the Tiller (U) to the lowest position. When the Tiller (U) is in the lowest position it can be locked by the pin of tiller angle adjustment.

Completion of disassembling:

The Scooter can be disassembled into 7 main parts.

Re-assembling your scooter by vice versa above disassembling procedures.

8. Caution: Obstacle climbing

1. Your scooter can climb obstacles and curbs of up to 50mm / 2" in height. Never attempt to overcome an obstacle when on an uphill or downhill gradient! Always approach obstacles straight on! Ensure that the front wheels and rear wheels move over the obstacle in one stroke, do not stop halfway!
2. The maximum gap the scooter can drive over is 100 mm / 4".

When driving scooter on ramp, adjust body center of gravity to keep scooter more safety.

Others:

1. Charge the batteries after each trip. If the scooter is not used for some time, batteries may lose their charge. Batteries should be charged at least once a month.
2. Check the battery gauge before driving to prevent power depletion.
3. Batteries will have an aging phenomenon, where the storage capacity will gradually decrease. If batteries are damaged, please wrap them in a plastic bag and contact your local dealer for proper disposal.
4. Do not disassemble battery and open sealed parts by yourself to prevent electric shock and burns from acid leakage.
5. Adjust speed to a slow speed when starting off to prevent sudden acceleration.
6. Never attempt to drive downhill backwards.
7. Try not to drive scooter at night or in rain or bad weather.
8. Your scooter is splash-proof. If the scooter has become wet from rain, dry the scooter before parking it for storage.
9. If storing your scooter for a long time (1 month or more), make ure that battery are fully charged, then disconnect the two batteries plugs (W), and store the scooter in a dry location.
10. Font basket, weight capacity 3kgs (6.5lbs)

9. Care and Maintenance

Daily Check

Please always check your scooter before you start your every journey.

GB

Check point	Inspection	What to do if the inspection is failed?
N-D lever	Check for correct function.	Contact your dealer.
Horn	Check for correct function.	Contact your dealer.
Throttles	Pull the wigwag to test the scooter movement.	Contact your dealer.
Electromechanical brake and Emergency hand brake	Pull the wigwag a little bit and release it to test if brake works. If your scooter comes with emergency handbrake, please check it as well.	Contact your dealer.
Battery Gauge	Check if the battery gauge is displayed and whether it is at low power.	<ul style="list-style-type: none"> • Contact your dealer if battery gauge is not working. • Recharge the battery immediately if low.
Rear mirror (s)	Check if the parts are clean and firmly tighten to the scooter and do not wobble.	<ul style="list-style-type: none"> • Clean up the dirt by damp cloth. • Tighten the screw or damping stem that holds the mirror(s)
Lighting	Check if all lights, such as head lights, rear lights, an turn signal are functioning correctly.	Contact your dealer.

Weekly Check

Check point	Inspection	What to do if the inspection is failed?
Speed Dial Knob	Check for correct function.	Contact your dealer.
Armrests	Check if the parts are clean and firmly tighten to the scooter and do not wobble. Tighten the screw knob that holds the armrest.	Contact your dealer.
Wheels/Tires	Inflate the tire to the correct pressure and check that: <ol style="list-style-type: none"> 1. Drive wheels rotate without wobbling. 2. Tire tread depth is not less than 0.5 mm. 3. No foreign objects in tires. 	Contact your dealer.
Motor	Check for correct function.	Contact your dealer.
Battery Charger	Check if the charger is functioning correctly and the batteries are charged.	Contact your dealer.

Monthly Check

Check point	Inspection	What to do if the inspection is failed?
Seat / Upholstery	Check for movement and if it's worn.	Contact your dealer.
Electronics	Check if all cables and connectors are firmly tighten to the scooter.	Contact your dealer.

Cleaning your Scooter:

- Do not use any abrasive or scouring liquids for cleaning. Only use a damp cloth and gentle detergent.
- Do not use hose pipe or splash water directly onto the scooter as this may cause damage to electronics.

GB

Maintenance:

- User should inspect the scooter regularly to keep scooter in good running order.
- Check if the electrical cable connectors are fully connected.
- All maintenance and repair of scooter should be done by an authorized dealer.

Seat Upholstery:

Only use damp cloth and a little soap to wipe the seat. Do not use abrasive cleaners as this will damage the seat.

Storage:

- Please store the scooter in a dry location. If store the scooter in long time, please disconnect the battery terminals.
- Do not store your scooter where it will be exposed to source of direct heat, damp, oil, acid, alkaline, or where Ozone could be possibly generated. All of the above will minimize scooter / tire cycle and shorten its lifetime.
- Store your scooter in a place where temperatures are between 10°C and 30°C.

10. Other Information:

Recycling & Disposal

- The equipment wrapping is potentially recyclable.
- The metal parts are used for scrap metal recycling. The plastic parts are used for plastic recycling.
- Electric components and printed circuit boards are disposed of as electronic scrap.
- Exhausted or damaged batteries can be returned to your medical equipment supplier.
- Disposal must be carried out in accordance with the respective national legal provisions.
- Ask your city or district council for details of the local waste management companies.

Service life:

We estimate a service life of five years for this product, provided it is used in strict accordance with the intended use as set out in this document and all maintenance and service requirements are met. The estimated service can be exceeded if the product is carefully used and properly maintained, and provided technical and scientific advances do not result in technical limitations. The service life can also be considerably reduced by extreme or incorrect usage. The fact that we estimate a service life for this product does not constitute an additional warranty.

Trouble Shooting:

Here are some suggestions about solving problems you may have with your scooter. There is a self-diagnostic warning light on the control panel. To check the self-diagnostic warning light, turn on the key and count the number of blinks on the warning light.

Scooter won't move when Key is turned on:

Check Point	Solution
Check if the power is off	Turn the ower on.
Check if the N-D lever is in Neutral position	Switch to D (drive) position. Turn off the power and turn on again.
Check if the battery power is enough. (Battery gauge is under 25%)	Recharge the battery and then retest.
Check if the charger power cord is still plugged in scooter	Unplug the charger power cord.

11. Error code:

GB

Flash	Description	Initial checkpoints
1	Battery low	<p>The batteries are running low.</p> <ul style="list-style-type: none"> Recharge the batteries.
2	Low Battery Fault	<p>The batteries have run out of charge.</p> <ul style="list-style-type: none"> Recharge the batteries.
3	High Battery Fault	<p>Battery voltage is too high. This may occur if overcharged &/or traveling down a long slope.</p> <ul style="list-style-type: none"> If traveling down a slope, reduce your speed to minimize the amount of regenerative charging. Check the battery and associated connections and wiring.
4	Current Limit Time-out	<p>The motor has been exceeding its maximum current rating for too long. This may be due to a faulty motor.</p> <ul style="list-style-type: none"> Check the motor and associated connections and wiring. Turn the controller off, leave for a few minutes and turn back on again.
5	Park Brake Fault	<p>Either a park brake release switch is active or the park brake is faulty.</p> <ul style="list-style-type: none"> Check the park brake and associated connections and wiring. Ensure any associated switches are in their correct positions.
6	Throttle	<p>The Throttle is out of neutral when turning the controller on.</p> <ul style="list-style-type: none"> Ensure the throttle is in neutral when turning the controller on. The Throttle may require re-calibration.

Flash	Description	Initial checkpoints
7	Speed Pot Fault	<p>The throttle, speed limit pot or their associated wiring may be faulty.</p> <ul style="list-style-type: none"> • Check the throttle and speed pot and associated connections and wiring.
8	Motor Voltage Fault	<p>The motor or its associated wiring is faulty.</p> <ul style="list-style-type: none"> • Check the motor and associated connections and wiring.
9	Other error	<p>The controller may have an internal fault.</p> <ul style="list-style-type: none"> • Check all the connections & wiring.

12. Specification

Specifications	REHASHOP electric scooter Siam
Overall Length	1020 mm / 40.2"
Overall Width	500 mm / 20"
Overall Height	860 mm / 34"
Front Wheels	200 mm / 8"
Rear Wheels	200 mm / 8"
Weight W/ Batteries	42,5 kg / 93.7 lbs
Weight of heaviest Piece	16 kg / 35.3 lbs
Max. Speed	6 km/h / 4 mp/h
Weight Capacity	136 kg / 300 lbs
Ground Clearance	40 mm
Grade Climbable	8 Grad
Curb Climbable	50 mm
Turning Radius	1060 mm
Brake	Electro-Mechanical
Seat Type	Swivel Padded Foldable
Seat Width	425 mm / 16.7"
Motor size	200W, 4650 r.p.m
Batterie size	(2) 12V. 12Ah
Weight of Battery Pack	9,5 kg / 20.9 lbs
Travel Range	11 km / 6.8 Miles
Battery Charger	2A (On Board)
Elektronics	On/off Key Switch, Battery Level Indicator, Speed Control Knob


Warranty Information:

The warranty period of the supplied batteries is six months.

The warranty period of the scooter is two years in total (except for the wearing parts).

13. Spare parts list REHASHOP electric scooter Siam:

Nr.	Art.-Nr.	Description
1	REHAET-0184	Basket for electric scooter Siam
2	REHAET-0185	Seat without armrests for electric scooter Siam
3	REHAET-0186	Armrest Right for electric scooter Siam
4	REHAET-0187	Armrest Left for electric scooter Siam
5	REHAET-0188	Battery Box (with Accu) for electric scooter Siam
6	REHAET-0189	Rear wheel with rim for electric scooter Siam
7	REHAET-0190	Charging Cable for electric scooter Siam
8	REHAET-0191	Front wheel with rim for electric scooter Siam

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PROTENO GmbH
Kaphingst Online GmbH